



All models in this series have the same basic features, lenses, magnification and methods of operation. The only differences are:

- 106, 107 & 108: Illumination is supplied by a plano-concave mirror.
- 106-L, 107-L, 108-L, 109-L: Illumination is supplied by a substage electric illuminator (110v 15 watt).
- 106, 106-L, 107 & 107-L: Have standard coarse focusing only.
- 108, 108-L, 109-L: Have both coarse and fine focusing.



Models 106, 107 & 108 supplied with substage mirror.



Models 108, 108-L & 109-L feature coarse and fine focusing.